

Jadex Japan

Operation Theater Lights

Medical Lamps

Operating Lamp, Camera System

Model Jadex-7KT



System enables to monitor the operation on the screen.

Lamp components are standart lamp head, camera lamp head and monitor.

Lamp Specifications:

- Diameter of Headgear : Ø695 mm
- Illumination : 155.000 Lux
- Color Temperature : 4300 K
- Depth of Illumination : 720 mm
- Focus Diameter : 120-300 mm
- Bulb : 24V / 150W halogen
- Inlet Voltage : 220-240 V/50-60 Hz

Camera Specifications:

- Image Sensor : 1/4" SONY HAD CCD
- Signal system : PAL
- Zoom : 10 X optic / 12 X digital
- Effective resolution : 440.000 pixel

Operating Lamp, Triple

Model Jadex-37T

- Diameter of Headgear : $\varnothing 695$ mm
- Illumination : 155.000 Lux
- Color Temperature : 4300 K
- Depth of Illumination : 720 mm
- Focus Diameter : 120-300 mm
- Bulb : 24V / 150W halogen
- Inlet Voltage : 220-240 V/50-60 Hz



Operating Lamp, Double

Model Jadex-77T

- Diameter of Headgear : $\varnothing 695$ mm
- Illumination : 155.000 Lux
- Color Temperature : 4300 K
- Depth of Illumination : 720 mm
- Focus Diameter : 120-300 mm
- Bulb : 24V / 150W halogen
- Inlet Voltage : 220-240 V/50-60 Hz





Operating Lamp, Single

Model Jadex-70T

- Diameter of Headgear : Ø695 mm
- Illumination : 155.000 Lux
- Color Temperature : 4300 K
- Depth of Illumination : 720 mm
- Focus Diameter : 120-300 mm
- Bulb : 24V / 150W halogen
- Inlet Voltage : 220-240 V/50-60 Hz



Operating Lamp, Mobile & Battery

Model Jadex-5M/Jadex-5MB

Specifications:

- Diameter of Headgear : Ø500 mm
- Illumination : 60.000 Lux
- Color Temperature : 4300 K
- Depth of Illumination : 700 mm
- Focus Diameter : 120 mm
- Bulb : 24V / 110W halogen
- Inlet Voltage : 220-240 V/50-60 Hz

- Battery system provides continuous energy for 2,5 hours
- Battery runs automatically in case of electricity failure
- LED indicator shows the battery status

Minor Surgery Lamp, Wall Type

ModelJadex-4D



Specifications:

- Diameter of Headgear : Ø380 mm
- Illumination : 60.000 Lux
- Color Temperature : 4300 K
- Depth of Illumination : 700 mm
- Bulb : 22V / 80W
- Inlet Voltage : 220-240 V/50-60 Hz

Minor Surgery Lamp, Ceiling Type

Model Jadex-4T



Specifications:

- Diameter of Headgear : Ø380 mm
- Illumination : 60.000 Lux
- Color Temperature : 4300 K
- Depth of Illumination : 700 mm
- Bulb : 22V / 80W
- Inlet Voltage : 220-240 V/50-60 Hz

Minor Surgery Lamp, Mobile & Battery

Model Jadex-4M / Jadex-4MB



Specifications:

- Diameter of Headgear : Ø380 mm
- Illumination : 60.000 Lux
- Color Temperature : 4300 K
- Depth of Illumination : 700 mm
- Bulb : 22V / 80W
- Inlet Voltage : 220-240 V/50-60 Hz

- Battery system provides continuous energy for 2 hours
- Battery runs automatically in case of electricity failure
- LED indicator shows the battery status

Magnifier Examination Lamp

Model Jadex-10



Specifications:

- Magnifier : 3,5 dioptries
- Magnifier size : 120 mm
- Bulb : 22 W circular fluorescent lamp
- Inlet voltage : 220-240V AC/50-60Hz



Examination Lamp, Mobile Model Jadex-40

Specifications:

- Illumination : 32.000 Lux / 0.5 m
- Lamp Head Dia. : 100 mm
- Bulb : 12V / 35W Halogen
- Inlet voltage : 220-240V AC/50-60Hz



Examination Lamp, Mobile with Flexible Arm

Model Jadex-60

Specifications:

- Illumination : 22.000 Lux / 0.5 m
- Lamp Head Dia. : 100 mm
- Bulb : 12V / 35W Halogen
- Inlet voltage : 220-240V AC/50-60Hz



Examination Lamp, Rail Type with Flexible Arm

Model Jadex-70

Specifications:

- Illumination : 22.000 Lux / 0.5 m
- Lamp Head Dia. : 100 mm
- Bulb : 12V / 35W Halogen
- Inlet voltage : 220-240V AC/50-60Hz